


FAQ: 0005: Material Analysis Report

1. Error: Cannot Continue: This ArNo is not latest certificate.



The screenshot shows the PharmaSuite interface for a Material Analysis Report. The main form contains the following data:

- Location: VATVA, Make: SOHAM, Material Type: RM
- Material Code: 1PARA01, PARACETAMOL
- Mat. Inw No.: DVDR090043, Mat. Inw Date: 24-07-2009, Stock Qty.: 20.0000 KG
- Test Slip No.: VDMTR90018, Date Of T.R.: 27-08-2009, Batch No.: 10, Sample Size: 0.0200 KG
- Supplier: SOHAM PHARMACEUTICALS PVT. LTD., Test as per: IP
- Mfg. Name: AAKASH MANTHAN INDUSTRIES, Mfg. Date: MAR-2009, ExpDt: FEB-2014
- Reject / Release Date: [Empty], A.R. No.: DVDR090043, Status: [Empty]
- Reason/Remark: [Empty], Analysis Date: [Empty]
- Assay: [Empty] %
- Analyst: [Empty], Partly Reject Qty: [Empty]
- Used for calculation while: [Empty], Retest Date: [Empty]

An error dialog box titled "Forms" is displayed in the center, with the message: "Cannot Continue : This ArNo is not latest certificate". The dialog has an "OK" button. A black arrow points from the error message to the "Create Report" button in the main form.

At the bottom of the form, there is a "Test" table with columns for "Test" and "Specifications/Limits". Below the table, there are fields for "Only Numeric Result" and "Comply" (set to "Yes").

The status bar at the bottom of the window displays the error message: "Cannot Continue : This ArNo is not latest certificate" and "Record: 1/?".

Reason for error: Q.C.Sample of respective inward no. is done two times.

Solution: Delete one record from

- Utility->Material Inw.No.Trekking, enter inward no. and double click on particular entry.